HBsAg

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are deviations from the instructions in this package insert.

NAME AND INTENDED USE

The Determine™ HBsAg is an *In Vitro*, visually read, qualitative immunoassay for the detection of Hepatitis B Surface Antigen (HBsAg) in human serum, plasma or whole blood. The test is intended as an aid to detect HBsAg from infected individuals.

SUMMARY AND EXPLANATION OF THE TEST

Assays for HBsAg are used to screen blood and blood products for the presence of HBsAg to prevent transmission of hepatitis B virus (HBV) to recipients of these products. HBsAg assays are also routinely used to diagnose suspected HBV infection and to monitor the status of infected individuals; i.e., whether the patient's infection has resolved or the patient has become a chronic carrier of the virus.¹

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

Determine HBsAg is an immunochromatographic test for the qualitative detection of Hepatitis B Surface Antigen (HBsAg).

Sample is added to the sample pad. As the sample migrates through the conjugate pad, it reconstitutes and mixes with the selenium colloid-antibody conjugate. This mixture continues to migrate through the solid phase to the immobilized antibodies at the patient window site.

If HBsAg is present in the sample, the antigen binds to the antibody-selenium colloid and to the antibody at the patient window, forming a red line at the patient window site.

If HBsAg is absent, the antibody-selenium colloid flows past the patient window, and no red line is formed at the patient window site.

To insure assay validity, a procedural control bar is incorporated in the assay device.

CONTENTS

Determine HBsAg Serum/Plasma Assay (List No. 7D25-43), 100 Tests

- Determine HBsAg Test Card, 10 cards (10 tests/card), Anti-HBs (mouse, monoclonal) antibody coated.

Determine HBsAg Whole Blood Assay (List No. 7D25-53), 100 Tests

- Determine HBsAg Test Card, 10 cards (10 tests/card), Anti-HBs (mouse, monoclonal) antibody coated.
- 1 Bottle (2.5 mL) Chase Buffer (List No. 7D22-43) prepared in phosphate buffer. Preservatives: Antimicrobial Agents.

ACCESSORIES (required but not provided)

Whole Blood (fingerstick assay)
EDTA Capillary Tubes No. 7D22-22

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use.

CAUTION:

Appropriate biosafety practices^{2,3} should be used when handling specimens and reagents. These precautions include, but are not limited to the following:

- Wear gloves.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of specimens or reagents using a suitable disinfectant, such as 0.5% sodium hypochlorite.^{4,5}
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local regulations.⁶

STORAGE

The Determine HBsAg Test Cards and Chase Buffer must be stored at 2-30°C until expiration date.

Kit components are stable until expiration date when handled and stored as directed. Do not use kit components beyond expiration date.

SPECIMEN COLLECTION

Serum, Plasma, and Whole Blood Collection by Venipuncture


Human serum, plasma, and whole blood collected by venipuncture should be collected aseptically in such a way as to avoid hemolysis.

NOTE: For whole blood specimens, EDTA collection tubes must be used.


Whole Blood Collection by Fingerstick⁷

Before collecting a fingerstick specimen, place an EDTA capillary tube on a clean dry surface.


1. Choose the fingertip of the middle, ring, or index finger (whichever is the least callused) for adults and children older than one year. Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.

2. Clean fingertip with alcohol; allow to air dry. Position the hand palm-side up.

3. Use a new lancet for each person. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose of the lancet in an appropriate biohazard sharps container.

4. Wipe away the first drop of blood with a sterile gauze pad.

5. Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times. Touch the tip of the EDTA Capillary Tube to the drop of blood*. Avoid air bubbles.

*If EDTA Capillary Tubes (No. 7D22-22) will be used, fill the tube with blood between the 2 marked lines.

SPECIMEN STORAGE

- Serum and plasma specimens should be stored at 2-8°C if the test is to be run within 7 days of collection. If testing is delayed more than 7 days, the specimen should be frozen (-20°C or colder).
- Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 7 days of collection. Do not freeze whole blood specimens.
- Whole blood collected by fingerstick should be tested immediately.

TEST PROCEDURE

The desired number of test units from the 10-test card can be removed by bending and tearing at the perforation.

NOTE: Removal of the test units should start from the right side of the test card to preserve the lot number which appears on the left side of the test card.

1. Remove the protective foil cover from each test.
2. For serum or plasma samples:
 - a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
 - b. Wait a minimum of 15 minutes (up to 24 hours) and read result.
3. For whole blood (venipuncture) samples:
 - a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
 - b. Wait one minute, then apply one drop of Chase Buffer to the sample pad.
 - c. Wait a minimum of 15 minutes (up to 24 hours) and read result.
4. For whole blood (fingerstick) samples:
 - a. Apply 50 µL of sample (by EDTA capillary tube) to the sample pad (marked by the arrow symbol).
 - b. Wait until blood is absorbed into the sample pad, then apply one drop of Chase Buffer to the sample pad.
 - c. Wait a minimum of 15 minutes (up to 24 hours) and read result.

QUALITY CONTROL

To insure assay validity, a procedural control is incorporated in the device and is labeled "Control". If the control bar does not turn red by assay completion, the test result is invalid and the sample should be retested.

INTERPRETATION OF RESULTS

POSITIVE (Two Bars)

Red bars appear in both the control window (labeled "Control") and the patient window (labeled "Patient") of the strip. Any visible red color in the patient window should be interpreted as positive.

NEGATIVE (One Bar)

One red bar appears in the control window of the strip (labeled "Control"), and no red bar appears in the patient window of the strip (labeled "Patient").

INVALID (No Bar)

If there is no red bar in the control window of the strip, and even if a red bar appears in the patient window of the strip, the result is invalid and should be repeated.

NOTES:

- The test result is positive even if the patient bar appears lighter or darker than the control bar.

LIMITATIONS OF THE PROCEDURE

- The Determine HBsAg test is designed to detect Hepatitis B Surface Antigen (HBsAg) in human serum, plasma, and whole blood. Other body fluids or pooled specimens may not give accurate results.
- The intensity of the patient bar does not necessarily correlate to the titer of antigen in the specimen.
- No test provides absolute assurance that a sample does not contain low levels of HBsAg such as those present at a very early stage of infection. Therefore, a negative result at any time does not preclude the possibility of exposure to or infection with HBsAg.
- Whole blood specimens containing anticoagulants other than EDTA may give incorrect results.
- For diagnostic purposes and to differentiate acute HBV infection from chronic HBV infection, the detection of HBsAg must be correlated with patient symptoms and other hepatitis B viral serological markers.

PERFORMANCE CHARACTERISTICS

SPECIFICITY

A total of 1,908 serum, plasma, and whole blood specimens from Asia, West Africa, and North America were tested by Determine HBsAg and Dainascreen HBsAg (Table I).

Population	Number of Specimens Tested	Negative by Determine HBsAg	Negative by Dainascreen HBsAg
Seronegatives			
Serum	682	681/682 (99.85%)	681/682 (99.85%)
Plasma*	498	498/498 (100.00%)	498/498 (100.00%)
Whole Blood	458	458/458 (100.00%)	458/458 (100.00%)
Pregnant Females	58**	57/57 (100.00%)	57/57 (100.00%)
West Africans	50	50/50 (100.00%)	50/50 (100.00%)
Disease States Other than HBV and Potentially Interfering Substances	162**	159/159 (100.00%)	159/159 (100.00%)
Total	1,908	1,903/1,904 (99.95%)	1,903/1,904 (99.95%)

*Ninety-two percent (458/498) of the plasma specimens were whole blood matched.

**One specimen from a pregnant female and three HIV-1 antibody positive specimens were positive by both Determine HBsAg and Dainascreen HBsAg. All four specimens confirmed positive by MEIA.

SENSITIVITY

A total of 434 HBsAg seropositive serum and whole blood specimens from Asia and North America were tested by Determine HBsAg and Dainascreen HBsAg (Table II).

Population	Number of Specimens Tested	Positive by Determine HBsAg	Positive by Dainascreen HBsAg
HBsAg Positive Serum	373	353/373 (94.64%)	336/373 (90.08%)
Whole Blood	61	60/61 (98.36%)	N/A N/A
Total	434	413/434 (95.16%)	336/373 (90.08%)

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BIBLIOGRAFIA (vide última página)

BIBLIOGRAPHY / LITERATUR / BIBLIOGRAPHIE / BIBLIOGRAFÍA / BIBLIOGRAFIA

1. Hollinger FB. Hepatitis B Virus. In: Fields BN, Knipe DM, et al., eds. *Virology*, 2nd edition. New York: Raven Press Ltd., 1990:2171-236.
2. World Health Organization. *Laboratory Biosafety Manual Third edition*. Geneva: World Health Organization, 2004.
3. CLINICAL AND LABORATORY STANDARDS INSTITUTE. *Protection of Laboratory Workers From Occupationally Acquired Infections*; Approved Guideline-Third Edition M29-A3 Vol. 25 No. 10
4. CDC. Recommendations for Prevention of HIV Transmission in Health-Care Settings, *MMWR* 1987; 36(SU02); 001
5. Sehulster LM, Hollinger FB, Dreesman GR, et al. Immunological and Biophysical Alteration of Hepatitis B Virus Antigens by Sodium Hypochlorite Disinfection. *Appl Environ Microbiol* 1981;42(5):762-7.
6. CLINICAL AND LABORATORY STANDARDS INSTITUTE. *Clinical Laboratory Waste Management*; Approved Guideline-Second Edition GP5-A2 Vol. 22 No. 3
7. CLINICAL AND LABORATORY STANDARDS INSTITUTE. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Fifth Edition H4-A5 Vol. 24 No. 21

Advice Line

For further information, please contact your distributor,
or call Inverness Medical Technical Specialists on:

UK: 08705 134952

International: +44 (0)1234 835959

<http://www.determinetest.com>

Mfd. for / Hergest. für / Fab. pour / Fab. para / Prod. per: Inverness Medical Japan Co., Ltd.
By / von / par / por / da: ABBOTT JAPAN CO., LTD.

March, 2007
©2007 Inverness Medical Japan Co., Ltd.