



- [Syphilis TP](#)
- [Features & Benefits](#)
- [Performance Data](#)
- [Determine Syphilis TP How to Use](#)
- [Determine Syphilis Product](#)
- [Specifications](#)
- [Determine Syphilis FAQ](#)



Features & Benefits

Determine Syphilis TP

Determine® Syphilis TP is an in vitro, visually read, qualitative immunoassay for the detection of antibodies to *Treponema pallidum* in human serum, plasma and whole blood.

With its simple one-step procedure for serum/plasma, two step procedure for whole blood, **Determine® Syphilis TP** is quick and easy to use, delivering clear, dependable results in just 15 minutes.

By diagnosing infection so rapidly, treatment can commence at the same visit, saving time and convenience for the health professional and the patient, as well as reducing the patient's risk of acquiring other sexually transmitted infections.

Scientifically engineered to provide dynamic performance in a rugged and robust format, **Determine® Syphilis TP** can be shipped, stored and used in nearly any environmental condition (2°C-30° C), even where there is no power or water.

Easy to use and highly flexible, **Determine® Syphilis TP** requires minimal training, no equipment and can be relied upon for fast, accurate diagnosis in the most basic of healthcare settings.

Refrigeration storage is not required with **Determine® Syphilis TP** and it has built in flexibility to be used with whole blood or serum/plasma, collected either by a fingerstick or venipuncture.

[Determine HIV Test](#) [Determine HIV-1/2 Ag/Ab](#) [Determine Hepatitis B Test](#) [Determine Syphilis Test](#) [Site Map](#)

Determine HIV-1/2 Ag/ Ab Combo is currently not available in the US