

- [HIV 1/2](#)
- [Features & Benefits](#)
- [Determine HIV 1/2 Performance Data](#)
- [Determine HIV 1/2 How to Use](#)
- [Determine HIV 1/2 Product](#)
- [Specifications](#)
- [Determine HIV 1/2 FAQ](#)

## Product Specifications



### Determine HIV 1/2

<b>Catalogue Number</b>	<b>Product</b>	<b>CE Marked</b>	
	Serum and Plasma 20 test	7D2342	7D2346
	Serum and Plasma 100 test	7D2343	7D2347
	* Whole Blood 20 test	7D2352	7D2346
	* Whole Blood 100 test	7D2353	7D2347
<b>Chase Buffer (CE and Non CE Marked)</b>		7D2243	
<b>EDTA capillaries</b>		7D2222 7D2227	
<b>Method</b>	Lateral flow		
<b>Time to result</b>	15 minutes		
<b>Storage conditions</b>	2-30°C		
<b>Shelf-life</b>	14 months		
<b>Sensitivity</b>	HIV-1 Positive 100% HIV-2 Positive 100% HIV-1 Subtype A-G 100% HIV-1 Group O 100% Europe 100% North America 100% Asia 100% Africa 99.91%		
<b>Specificity</b>	99.75%		
<b>Sample Type</b>	Serum/plasma/whole blood		
<b>Kit size</b>	20 and 100 tests		
<b>Kit contents</b>	20 or 100 test units		

\* For Whole Blood please order Chase Buffer

[Determine HIV Test](#) | 
 [Determine HIV-1/2 Ag/Ab](#) | 
 [Determine Hepatitis B Test](#) | 
 [Determine Syphilis Test](#) | 
 [Site Map](#)

Determine HIV-1/2 Ag/ Ab Combo is currently not available in the US

