



- HBsAg
- Features & Benefits
- Determine HBsAg Performance Data
- How to Use
- Determine HBsAg Product
- Specifications
- Determine HBsAg FAQ



Product Specifications



Determine HBsAg

Catalogue Number	Product	
	Serum and Plasma 100 test	7D2543
	* Whole Blood 100 test	7D2553
Chase Buffer		7D2243
EDTA capillaries		7D2222
Method	Lateral flow	
Time to result	15 minutes	
Storage conditions	2-30°C	
Shelf-life	18 months	
Sensitivity	95.16%	
Specificity	99.95%	
Sample Type	Serum/plasma/whole blood	
Kit size	100 tests	
Kit contents	100 test units	

** For Whole Blood please order Chase Buffer*

[Determine HIV Test](#) |
 [Determine HIV-1/2 Ag/Ab](#) |
 [Determine Hepatitis B Test](#) |
 [Determine Syphilis Test](#) |
 [Site Map](#)

Determine HIV-1/2 Ag/ Ab Combo is currently not available in the US