

- [Syphilis TP](#)
- [Features & Benefits](#)
- [Performance Data](#)
- [Determine Syphilis TP How to Use](#)
- [Determine Syphilis Product](#)
- [Specifications](#)
- [Determine Syphilis FAQ](#)

Product Specifications

Determine Syphilis TP



Catalogue Number	Product	CE and non-CE Marked	CE Marked
	Serum and Plasma 100 test	7D2443	
	* Whole Blood 100 test	7D2453	
Chase Buffer (CE and Non CE Marked)		7D2243	
EDTA capillaries		7D2222	7D2227
Method	Lateral flow		
Time to result	15 minutes		
Storage conditions	2-30°C		
Shelf-life	12 months		
Sensitivity	Whole Blood 92.31% Plasma 100% Serum 100%		
Specificity	Whole Blood 100% Plasma 100% Serum 100%		
Sample Type	Serum/plasma/whole blood		
Kit size	100 tests		
Kit contents	100 test units		

* For Whole Blood please order Chase Buffer

[Determine HIV Test](#) |
 [Determine HIV-1/2 Ag/Ab](#) |
 [Determine Hepatitis B Test](#) |
 [Determine Syphilis Test](#) |
 [Site Map](#)

Determine HIV-1/2 Ag/ Ab Combo is currently not available in the US