



- [About HIV](#)
- [Features & Benefits](#)
- [Determine Performance Data](#)
- [How to Use](#)
- [Determine Product Specifications](#)
- [Determine FAQ](#)



Features & Benefits

Determine® HIV-1/2 Ag/ Ab Combo

Determine® HIV-1/2 Ag/Ab Combo is a revolutionary, fourth generation, rapid *in vitro* immunoassay for the qualitative simultaneous detection of HIV p24 antigen as well as antibodies to HIV-1 and HIV-2 in serum, plasma and whole blood.

The p24 antigen (Ag) is produced during the first few weeks of HIV infection and is detectable before HIV antibodies (Ab) are produced, making it an ideal marker in aiding early HIV diagnosis. Due to separate signals for HIV Ag and anti-HIV Ab, **Determine® HIV-1/2 Ag/Ab Combo** provides determination and differentiation of the infection phase and enables the identification of acute HIV infection. Based on testing standard seroconversion panels, **Determine® HIV-1/2 Ag/Ab Combo** test detected HIV infection several days earlier than the HIV-antibody only product. **Determine® HIV-1/2 Ag/Ab Combo** is quick and easy to use, delivering clear, dependable results in just 20 minutes.

By detecting both p24 antigen and antibodies to HIV, **Determine® HIV-1/2 Ag/Ab Combo** will help healthcare workers around the globe identify individual infection, ideally preventing further transmission from those who may not be aware of their HIV status.

Scientifically engineered to provide dynamic performance in a rugged and robust format, **Determine® HIV-1/2 Ag/Ab Combo** can be shipped, stored and used in broad and challenging environmental conditions (2°C -30°C), even where there is no power or water.

Easy to use and highly flexible, **Determine® HIV-1/2 Ag/Ab Combo** requires minimal training, no equipment and can be relied upon for fast, accurate diagnosis in the most basic of healthcare settings.

Refrigeration storage is not required with **Determine® HIV-1/2 Ag/Ab Combo** and it has the built in flexibility to be used with whole blood or serum/plasma, collected either by a fingerstick or venipuncture.

[Determine HIV Test](#) [Determine HIV-1/2 Ag/Ab](#) [Determine Hepatitis B Test](#) [Determine Syphilis Test](#) [Site Map](#)

Determine HIV-1/2 Ag/ Ab Combo is currently not available in the US