

- [HBsAg](#)
- [Features & Benefits](#)
- [Determine HBsAg Performance Data](#)
- [How to Use](#)
- [Determine HBsAg Product](#)
- [Specifications](#)
- [Determine HBsAg FAQ](#)



Features & Benefits

Determine HBsAg

Determine® HBsAg is an in vitro, visually read, qualitative immunoassay for the detection of Hepatitis B Surface Antigens in human serum, plasma and whole blood.

With its simple one-step procedure for serum/plasma, two step procedure for whole blood, **Determine® HBsAg** is quick and easy to use, delivering clear, dependable results in just 15 minutes.

By detecting early markers of infection, **Determine® HBsAg** will help diagnose acute infection, confirm chronic infection and screen blood supplies allowing healthcare professionals to initiate prompt and appropriate follow-up action to help infected patients.

Scientifically engineered to provide dynamic performance in a rugged and robust format, **Determine® HBsAg** can be shipped, stored and used in nearly any environmental condition (2°C-30° C), even where there is no power or water.

Easy to use and highly flexible, **Determine® HBsAg** requires minimal training, no equipment and can be relied upon for fast, accurate diagnosis in the most basic of healthcare settings.

Refrigeration storage is not required with Determine® HBsAg and it has built in flexibility to be used with whole blood or serum/plasma, collected either by a fingerstick or venipuncture.

[Determine HIV Test](#) [Determine HIV-1/2 Ag/Ab](#) [Determine Hepatitis B Test](#) [Determine Syphilis Test](#) [Site Map](#)

Determine HIV-1/2 Ag/ Ab Combo is currently not available in the US